

How to add/setup a new tag

Ivana Ivancic - 2021-03-10 - Blyott 'How To' Guides



This page will provide step by step instructions on how to add/setup a new tag attached to a device for tracking within the Smart Healthcare platform.

Instructions

1. Select the 'Tags' tab from the left-hand panel (take note that different screen sizes change the layout of the panel, on smaller screens the panel will be on the right side because of the responsiveness of the webpage - RWD).
2. Select the '+' icon from the top right corner of the screen.
3. Complete all data fields.
4. Be sure to enter the correct MAC ID/NFC ID. The number will not be visible on the tag, but can easily be checked using a mobile app downloaded from your relevant app store.
5. Be sure to use the appropriate name and description.
6. Be sure to use the correct (local) asset code.

Once tags have been added, they can be searched for within the platform by selecting the 'filter' icon from the top right corner of the screen.

Only a user with 'Administration' access to the platform will be able to add/setup a new Tag.

If you have any questions or suggestions at all, please feel free to get in touch with us [HERE!](#)

Thanks and regards,

The Blyott Support Team



www.blyott.com | support@blyott.com | support.blyott.com